



Fraunhofer

TESTED[®] DEVICE

Memmert GmbH + Co. KG
Heating cabinet UF 750 plus
Report No. ME 1302-636

Qualification Certificate

This is to certify that the product mentioned above, provided by

Memmert GmbH + Co. KG
Schwabach, Germany

has been awarded a Fraunhofer certificate TESTED DEVICE bearing the report number ME 1302-636.

Operated at 125 °C with the fan running at full speed, the universal heating cabinet UF 750 plus is suitable for use in cleanrooms fulfilling the specifications of Class 7 in accordance with ISO 14644-1. If the heating cabinet is supplied with Class 1 air, the quality of the air inside the cabinet corresponds with Class 7 in accordance with ISO 14644-1.

Operated at 125 °C with the fan at rest, the universal heating cabinet UF 750 plus is suitable for use in cleanrooms fulfilling the specifications of Class 5 in accordance with ISO 14644-1. If the heating cabinet is supplied with Class 1 air, the quality of the air inside the cabinet corresponds with Class 3 in accordance with ISO 14644-1.

Detailed information and parameters of the test environment can be found in the Fraunhofer IPA test report.

Stuttgart, February 14, 2013

i.A. 
Project manager

DUPLICATE

The validity of this certificate applies only to the mentioned product in this particular condition for a duration of 5 years.

Further information:
www.tested-device.com